




Product	Mattress
Item number	R8487378
Color	Navy blue
Product description	<p>VVS: 777541708 EAN: 5708590360560 HMI: 125838</p>  <p>Mattress for change table, 31.7" x 71.7". Incl. screws. The screws ensure that the mattress is firmly attached to the changing table – both during transfer to and from the changing table and when the changing table is folded away.</p> <p>Mattress for change table, 31.7" x 71.7". Incl. screws. The screws ensure that the mattress is firmly attached to the changing table – both during transfer to and from the changing table and when the changing table is folded away.</p>
Materials	The mattress cover is manufactured from 100% polyester with polyurethane, approved in accordance with Oeko-Tex 100, Class I. It is fire retardant and impermeable to dust mites. The foam filling is flexible polyester-based polyurethane foam, which is also fire retardant and approved in accordance with Oeko-tex 100, Class I.
Ergonomy	The mattress scan be turned over and is made from a pleasant hypoallergenic material with a surface that is easy to wipe down and keep clean.
Features	The mattress is easily removed, which offers flexibility of use and ensures optimum opportunity for effective cleaning over and under the mattress..



Product	Mattress
Item number	R8487378
Cleaning	<p>This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged.</p> <p>The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container.</p> <p>DISINFECTANTS: Sprit (husholdningssprit) Tvættsprit 70% Safe Zone Plus Selgiene Extreme Accel PREvention Wipes Lingettes OxiGenesis from ENVIROX</p> <p>The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.</p>