




# PRESSALIT®

|                                      |                                                                                                                                                                                                                                                                                                                                      |
|--------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product</b><br><b>Item number</b> | Mattress for MCT 2, MCT 3 and MCT 5<br>R842719378                                                                                                                                                                                                                                                                                    |
| <b>Colour</b>                        | VVS: 777542738<br>EAN: 5708590389837<br>HMI: 139823<br><br>Navy blue                                                                                                                                                                                                                                                                 |
| <b>Product description</b>           | <br><br>Mattress for MCT 2, MCT 3 and MCT 5 change trolleys, 695 x 1880 mm.                                                                                                                                                                       |
| <b>Materials</b>                     | The mattress cover is manufactured from 100% polyester with polyurethane, approved in accordance with Oeko-Tex 100, Class I. It is fire retardant and impermeable to dust mites. The foam filling is flexible polyester-based polyurethane foam, which is also fire retardant and approved in accordance with Oeko-tex 100, Class I. |
| <b>Warranty</b>                      | Pressalit makes a commitment to remedy material and manufacturing defects for a period of three years from the invoicing date.                                                                                                                                                                                                       |



# PRESSALIT®

|                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product</b>     | Mattress for MCT 2, MCT 3 and MCT 5                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Item number</b> | R842719378                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Cleaning</b>    | <p>CT and MCT mattresses</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged.</p> <p>The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container.</p> <p>DISINFECTANTS:</p> <p>Sprit (husholdningssprit)</p> <p>Tvättssprit 70%</p> <p>Safe Zone Plus</p> <p>Selgiene Extreme</p> <p>Accel PREvention Wipes Lingettes</p> <p>OxiGenesis from ENVIROX</p> <p>The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.</p> |

[pressalit.com](http://pressalit.com)

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

## PRESSALIT®