




| | |
|---|--|
| <p>Product</p> <p>Item number</p> | <p>TMA 1 arm support for mounting on toilet</p> <p>R1171000</p> |
| <p>Color</p> <p>Product description</p> | <p>VVS:615723100</p> <p>NRF:6191651</p> <p>EAN:5708590211718</p> <p>HMI: 42230</p> <p>NBS:N13:436</p> <p>DB: 2076624</p> <p>White</p>  <p>TMA arm supports comes in various widths, and models with or without toilet seat.</p> <p>TMA 1 arm supports for mounting on toilet, incl. BL1999 hinge. Dark gray handles. BL1 hinge, is only recommended with toilet seats R26, R27, R28, R29, R37 or R38.</p> |
| <p>Adjustability</p> | <p>Fixed mounting</p> |
| <p>Materials</p> <p>Ergonomy</p> <p>Features</p> | <p>Powder-coated stainless steel.</p> <p>Soft polyurethane.</p> <p><i>Flammability-tested (ISO 8191-1). Test report available.</i></p> <p>The folding arms provide firm support and the molded handles are comfortable to hold on to.</p> <p>The arm supports can be folded individually.</p> |
| <p>Warranty</p> <p>Cleaning</p> <p>Declaration of conformity</p> | <p>Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.</p> <p>This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged.</p> <p>The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.</p> <p>This product is in compliance with all relevant requirements of regulation (EU) 2017/745 on medical devices.</p> |



PRESSALIT®

| | |
|--------------------|--|
| Product | TMA 1 arm support for mounting on toilet |
| Item number | R1171000 |

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®