




PRESSALIT®

| | |
|----------------------------|---|
| Product | Shower mattress for MSCT 1 and MSCT 4 |
| Item number | R844721304 |
| | VVS: 777542821 EAN: 5708590373935 HMI: 128434 |
| Colour | Pearl grey |
| Product description | <p>Shower mattress for MSCT 1 and MSCT 4 shower change trolley, 660 x 2060 mm. Incl. 1.6 m hose.</p> <p>Shower mattress for MSCT 1 and MSCT 4 shower change trolley, 660 x 2060 mm. Incl. 1.6 m hose.</p> |
| Materials | The mattress cover is welded polyester fabric with PVC-coating. It is UV-resistant and 100% impermeable to water. It contains neither CFC nor latex. The foam filling is flexible polyether-based polyurethane foam, approved in accordance with Oeko-tex 100, Class I. |
| Warranty | Pressalit makes a commitment to remedy material and manufacturing defects for a period of three years from the invoicing date. |





| | |
|----------------------------------|---|
| Product | Shower mattress for MSCT 1 and MSCT 4 |
| Item number | R844721304 |
| Cleaning | <p>MSCT mattresses The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged.</p> <p>The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container.</p> <p>DISINFECTANTS: Sprit (husholdningssprit) Tvättspirit 70% Safe Zone Plus Selgiene Extreme Microsan Accel PREvention Wipes Lingettes</p> <p>The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.</p> |
| Declaration of conformity | <p>MDR-ACC This product fulfils all relevant requirements in Regulation (EU) 2017/745 on medical devices in risk class I. This product is considered to be an accessory for a medical device, as it is not itself a medical device, but is intended to be used together with several particular medical devices to specifically enable the medical devices to be used in accordance with their intended purposes.</p> |