

PRESSALIT.

Product	SCT 2000 shower change table, mesh cover, electrically height adjustable		
Item number	R8591218000		
	EAN:5708590376127		
Color	FDA: INQ/D124319 White		
	34.65 35.83 70.87		
Product description	Pressalit offers a variety of mobile and wall-hung, height-adjustable and fixed he shower and changing tables of high quality. They are designed to provide a stab platform for showering and changing adults and children with reduced mobility. These tables can be used in schools, institutions, hospitals and private homes a are all designed to be durable, comfortable and easy to clean. The design and flexibility of the tables ensures optimal working conditions for the carer as well as a high level of comfort and safety for the user. IAPMO-certified. SCT 2000 shower change table, powered, height adjustable 27.6"; from 11.8" - 3 with wired hand control Foldable. Incl. safety rail and water collection tray with integrated outlet. Outer measurements 35.8" x 70.9". Excl. mounting kit.		
Adjustability	Height adjustable		
Adjustment height	27,58 inch		
Product capacity	Maximum user weight (MUM): 441 lbs. Maximum load (SWL): 441 lbs. When the table is loaded with 441 lbs, each mounting bolt supports 57 lbs. This product has been tested statically with a load of 1102 lbs as the IEC 60601-1 standard requires testing with a safety factor 2.5.		
Speed up	Speed up: 0,75 inch/s Speed down: 0,83 inch/s		
Noise level	50 dB		



PRESSALIT.®

Product	SCT 2000 shower change table, mesh cover, electrically height adjustable	
Item number	R8591218000	
Materials	Construction in powder-laquered stainless steel, mesh cover in PVC-coated po (for extra durability), water collection tray in PVC-U and carriage housing in polystyrene.	
Ergonomy	The 2000 series changing table has a low start position to make transfers to the tal easier. The washable, perforated mesh cover lets the water pass through to the was collection tray underneath. The flexibility of the changing table provides for optimal working conditions for the carer as well as a high level of comfort and safety for the user.	
Features	A water collection tray with integrated water outlet prevents water from splashing or caregiver or the floor. The flexible hose attached to the center mounted outlet of the water collection tray, can be connected to either an outlet in the floor or on the wall. The integrated pneumatic cylinders enable the table to be folded out of the way who not in use to maximize floor space in front of the unit. Electric with wired hand control – the red buttons clearly indicate "Up" and "Down". The wired hand control comes with a holder for fitting to the wall. The stretcher of the table is very easy to clean, and the stretcher cover can be unzipped for washing or easy replacement.	
Mains connection (V)	120	
Actuator (motor)	LA40 Careline (Linak), 4000 N, IPX6	
Control box	CA63 Careline, (Linak), IPX6	
Wired remote control	HB30 (Linak), IP66	
Power consumption (maximum load)	Power consumption by maximum load: 420 W Power consumption by standby: 0,5 W	
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of three years from the invoicing date. NOTICE If the product is placed in an atmosphere heavy in salt or chloride, the warranty will be reduced from the normal 3 years to 1 year. Salt and chloride have a corrosive effect on the products components and will cause a decrease in the life expectancy of the product.	



PRESSALIT.

Product	SCT 2000 shower change table, mesh cover, electrically height adjustable		
ltem number	R8591218000		
Cleaning	This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged.		
	The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container.		
	DISINFECTANTS: Clinimax Difficil-S Erisan OXY+ Sprit (husholdningssprit) Suma Bac D10 Virkon Rodalon Desinfektionsspray		
	Lysoformin Spezial CaviCide Sufanios Premium Tvättsprit 70% Accel PREVention Wipes Lingettes OxiGenesis from ENVIROX		
	NON-DISINFECTANTS: Suma Star free D1 Taski 100 Taski 200 Milizid Citro		
	The product must be wiped down with a soft rag wrung out in clear water immedia after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.		



PRESSALIT®

Product		
Item number		

SCT 2000 shower change table, mesh cover, electrically height adjustable R8591218000

Mounting in different wall types

Mounting in different wall types			
	Concrete >B25 V8608 - Mounting kit, 6 lag bolts Ø10 x 70 - plugs SX12. To be ordered separately.		
	Solid Brick >= Mz12 V8608 - Mounting kit, 6 lag bolts Ø10 x 70 - plugs SX12. To be ordered separately.		
	Solid Sand Lime Brick >= HIz12 V8608 - Mounting kit, 6 lag bolts Ø10 x 70 - plugs SX12. To be ordered separately.		
	Aerated Concrete P4 V8658 - Mounting kit, 6 lag bolts Ø10 x 70 - plugs GB14. To be ordered separately.		
	Drywall (2 x 12,5 mm) water-resistant, re-enforced with wood V8608 - Mounting kit, 6 lag bolts Ø10 x 70 - plugs SX12. NB! Plugs are not needed. To be ordered separately.		



PRESSALIT®

Product Item number	SCT 2000 shower change table, mesh cover, electrically height adjustable R8591218000
	Drywall (2 x 12,5 mm) water-resistant, re-enforced with 22 mm plywood Wall fasteners (6 pcs) must be bought separately, for drywall it is possible to use e.g. SnapTight or WingIt anchors or similar.
	Drywall (2 x 12,5 mm) water-resistant, re-enforced with steel plates (.08 inch) Wall fasteners (6 pcs) must be bought separately, for drywall it is possible to use e.g. SnapTight or WingIt anchors or similar.
	Drywall (2 x 12,5 mm) water-resistant
	Mounting in this wall type is not recommended - risk of fatal injury to persons and walls.

It is the responsibility of the installer to ensure that suitable wall fasteners are used.

pressalit.com