




Product	Mattress
Item number	R8487378
	VVS: 777541708 EAN: 5708590360560 HMI: 125838
Colour	Navy blue
Product description	Mattress for change table, 805 x 1820 mm. Incl. screws. The screws ensure that the mattress is firmly attached to the change table – both during transfer to and from the changing table and when the changing table is folded away.  Mattress for change table, 805 x 1820 mm. Incl. screws. The screws ensure that the mattress is firmly attached to the changing table – both during transfer to and from the changing table and when the changing table is folded away.
Materials	The mattress cover is manufactured from 100% polyester with polyurethane, approved in accordance with Oeko-Tex 100, Class I. It is fire retardant and impermeable to dust mites. The foam filling is flexible polyester-based polyurethane foam, which is also fire retardant and approved in accordance with Oeko-tex 100, Class I. The mattress cover is manufactured from 100% polyester with polyurethane, approved in accordance with Oeko-Tex 100, Class I. It is fire retardant and impermeable to dust mites. The foam filling is flexible polyester-based polyurethane foam, which is also fire retardant and approved in accordance with Oeko-tex 100, Class I.
Ergonomy	The mattress scan be turned over and is made from a pleasant hypoallergenic material with a surface that is easy to wipe down and keep clean.
Features	The mattress is easily removed, which offers flexibility of use and ensures optimum opportunity for effective cleaning over and under the mattress.



PRESSALIT®

Product	Mattress
Item number	R8487378
Cleaning	CT and MCT mattresses The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged. The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container. DISINFECTANTS: Sprit (husholdningssprit) Tvättssprit 70% Safe Zone Plus Selgiene Extreme Accel PREvention Wipes Lingettes OxiGenesis from ENVIROX The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®