




<b>Product</b>	Mattress
<b>Item number</b>	R8487378
	VVS: 777541708 EAN: 5708590360560 HMI: 125838
<b>Colour</b>	Navy blue
<b>Product description</b>	<p>Mattress for change table, 805 x 1820 mm. Incl. screws. The screws ensure that the mattress is firmly attached to the change table – both during transfer to and from the changing table and when the changing table is folded away.</p> <p>Mattress for change table, 805 x 1820 mm. Incl. screws. The screws ensure that the mattress is firmly attached to the changing table – both during transfer to and from the changing table and when the changing table is folded away.</p>
<b>Materials</b>	<p>The mattress cover is manufactured from 100% polyester with polyurethane, approved in accordance with Oeko-Tex 100, Class I. It is fire retardant and impermeable to dust mites. The foam filling is flexible polyester-based polyurethane foam, which is also fire retardant and approved in accordance with Oeko-tex 100, Class I. The mattress cover is manufactured from 100% polyester with polyurethane, approved in accordance with Oeko-Tex 100, Class I. It is fire retardant and impermeable to dust mites. The foam filling is flexible polyester-based polyurethane foam, which is also fire retardant and approved in accordance with Oeko-tex 100, Class I.</p>
<b>Ergonomy</b>	The mattress scan be turned over and is made from a pleasant hypoallergenic material with a surface that is easy to wipe down and keep clean.
<b>Features</b>	The mattress is easily removed, which offers flexibility of use and ensures optimum opportunity for effective cleaning over and under the mattress.





# PRESSALIT®

<b>Product</b>	Mattress
<b>Item number</b>	R8487378
<b>Cleaning</b>	<p>CT and MCT mattresses</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged.</p> <p>The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container.</p> <p>DISINFECTANTS:</p> <p>Sprit (husholdningssprit)</p> <p>Tvåttsprit 70%</p> <p>Safe Zone Plus</p> <p>Selgiene Extreme</p> <p>Accel PREvention Wipes Lingettes</p> <p>OxiGenesis from ENVIROX</p> <p>The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.</p>

[pressalit.com](http://pressalit.com)

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

## PRESSALIT®