



# PRESSALIT®

|                                      |                                                                                                                                                                                                                                                                                                                                                                                                       |
|--------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product</b><br><b>Item number</b> | TMA 1 arm support for mounting on toilet<br>R1170000                                                                                                                                                                                                                                                                                                                                                  |
| <b>Color</b>                         | VVS:615723200<br>NRF:6191649<br>RSK:7828180<br>EAN:5708590211701<br>HMI: 35270<br>NBS:N13:436<br>FDA:ILS/D135866<br>White                                                                                                                                                                                                                                                                             |
| <b>Product description</b>           | TMA arm supports comes in various widths, and models with or without toilet seat.<br><br>TMA 1 arm supports for mounting on toilet, incl. R37 Dania toilet seat.<br>With BL1999 hinge. Dark gray handles.                                                                                                                                                                                             |
| <b>Adjustability</b>                 | Fixed mounting                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Materials</b>                     | Powder-coated stainless steel.<br>Soft polyurethane.<br><i>Flammability-tested (ISO 8191-1). Test report available.</i><br>Toilet seat and cover in thermoset plastic (urea).                                                                                                                                                                                                                         |
| <b>Ergonomy</b>                      | The folding arms provide firm support and the molded handles are comfortable to hold on to.                                                                                                                                                                                                                                                                                                           |
| <b>Features</b>                      | The arm supports can be folded individually.                                                                                                                                                                                                                                                                                                                                                          |
| <b>Warranty</b>                      | Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.                                                                                                                                                                                                                                                                         |
| <b>Cleaning</b>                      | This product should be cleaned with hot water or a mild cleaning solution.<br>Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged.<br><br>The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate. |



# PRESSALIT®

|                    |                                          |
|--------------------|------------------------------------------|
| <b>Product</b>     | TMA 1 arm support for mounting on toilet |
| <b>Item number</b> | R1170000                                 |

---

[pressalit.com](http://pressalit.com)

---

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

# PRESSALIT®