




# PRESSALIT®

<b>Product</b> <b>Item number</b>	TMA 3 arm support for mounting on toilet, with nurse call button to the right R1193112
<b>Colour</b>	Anthracite grey
<b>Product description</b>	<p>EAN:5708590387482</p>  <p>TMA arm supports comes in various widths, and models with or without toilet seat.</p> <p>TMA 3 arm supports for mounting on toilet. Anthracite grey handles. Incl. two paper holders. Excl. hinges. Bottom mounting. Compatible with IFØ Spira toilet 6261, Spira 6293 and Sign 6893 Recommended with toiletseats: Pressalit Spira V 1136011-DN3999 Pressalit Spira V 1136112-DN3999 Pressalit Sign 553011-DN1999 Pressalit Sign 553112-DN1999 Pressalit Sign 554011-DN1999 Pressalit Sign 554108-DN1999 Pressalit Sign 554112-DN1999</p>
<b>Adjustability</b> <b>Product capacity</b>	<p>Fixed mounting</p> <p>Maximum user weight (MUM): 150 kg. The armsupports are tested according to ISO/FDIS 17966:2025 with a safety factor of 1.5. The armsupports are additionally statically tested individually with 97 kg.</p>
<b>Materials</b> <b>Ergonomy</b> <b>Features</b> <b>Mains connection (V)</b>	<p>Anodised aluminium and ABS</p> <p>The arm supports provide firm support and the moulded handles are comfortable for the hands and arms.</p> <p>The arm supports are foldable. The arm supports can be folded up and down individually.</p> <p>40</p>



# PRESSALIT®

<b>Product</b>	TMA 3 arm support for mounting on toilet, with nurse call button to the right
<b>Item number</b>	R1193112
<b>Power consumption (maximum load)</b>	Power consumption by maximum load, watt: 4 W
<b>Warranty</b>	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.
<b>Cleaning</b>	<p><b>GENERAL</b></p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged.</p> <p>The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.</p>

[pressalit.com](http://pressalit.com)

# PRESSALIT®