




<b>Product</b> <b>Item number</b>	TMA 1 arm support for mounting on toilet R1171000
<b>Colour</b>  <b>Product description</b>	<p>VVS:615723100          NRF:6191651          EAN:5708590211718          HMI: 42230          NBS:N13:436          DB: 2076624</p> <p>White</p>  <p>TMA arm supports comes in various widths, and models with or without toilet seat.</p> <p>TMA 1 arm supports for mounting on toilet, incl. BL1999 hinge. Dark grey handles. BL1 hinge, is only recommended with toilet seats R26, R27, R28, R29, R30, R36, R37 or R38.</p>
<b>Adjustability</b>	Fixed mounting
<b>Materials</b>	Powder-coated stainless steel. Soft polyurethane. <i>Flammability-tested (ISO 8191-1). Test report available.</i>
<b>Ergonomy</b>	The folding arms give a good support and the moulded handles are comfortable for the hands and arms.
<b>Features</b>	The arm supports can be folded individually.
<b>Warranty</b>	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.
<b>Cleaning</b>	<p>General.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged.</p> <p>The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.</p>



PRESSALIT®

<b>Product</b>	TMA 1 arm support for mounting on toilet
<b>Item number</b>	R1171000
<b>Declaration of conformity</b>	CE UK MDR This product is in compliance with all relevant requirements of: - Regulation (EU) 2017/745 on medical devices. <b>For the UK:</b> - Medical Device Regulation (MDR) 2002 on medical devices.

---

[pressalit.com](http://pressalit.com)

---

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

PRESSALIT®