

PRESSALIT.

| Product | VALUE folding shower seat, fixed height |
|---------------------|---|
| Item number | R1600000 |
| | VVS:777526335 |
| | NRF:6191656 |
| | EAN:5708590176734 |
| | HMI: 25389 |
| | NBS:N13:448 |
| Colour | White 440 |
| Product description | VALUE folding shower seat and backrest create a simple solution for users who need to sit down while showering. |
| | Folding shower seat, seat pale grey Incl. mounting kit for brick wall |
| Adjustability | Fixed mounting |
| Product capacity | Max. vertical load: 125 kg. |
| Materials | Diecast, powder-coated aluminium. Rilsan-coated steel, which is comfortable to the touch, easy to clean, and protects the metal against corrosion. Plastic part is made of vacuum-moulded polyethylene (PE). Flammability-tested (ISO 8191-1 and ISO 8191-2). Test report available. Toilet seat and cover in thermoset plastic (urea). |
| Ergonomy | Ergonomic seat with a comfortable surface. |
| Features | The seat can fold up when not in use or to make more space. |
| | With run-off grooves. |
| Warranty | Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date. |



PRESSALIT.®

| Product | VALUE folding shower seat, fixed height |
|---------------------------|---|
| Item number | R1600000 |
| Cleaning | |
| | VALUE shower seat. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use corrosive or abrasive cleaners as it may result in the surface being damaged. For surface disinfection, chlorine-based disinfectants may be used according to the product specification and safety precautions for these products. |
| | The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container. |
| | DISINFECTANTS: Clinimax Difficil-S Erisan OXY+ Sprit (husholdningssprit) Suma Bac D10 Virkon Rodalon Desinfektionsspray Lysoformin Spezial CaviCide Sufanios Premium Tvättsprit 70% WipeClean UniversalCleaning Chrolorine Disinfection Wipes Klorilli (Kiilto Clean) |
| | NON-DISINFECTANTS: Suma Star free D1 Taski 100 Taski 200 Milizid Citro |
| | The product must be wiped down with a soft rag wrung out in clear water immediate after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate. |
| Declaration of conformity | CE MDR This product is in compliance with all relevant requirements of: - Regulation (EU) 2017/745 on medical devices. |



PRESSALIT®

| Duadriat | VALUE folding shower cost fixed height |
|-------------|---|
| Product | VALUE folding shower seat, fixed height |
| Item number | R1600000 |

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com