




# PRESSALIT®

<b>Product</b> <b>Item number</b>	VALUE backrest, fixed height R1651000
<b>Colour</b>	VVS:777526390 NRF:6191661 EAN:5708590176802 HMI: 25574 NBS:N13:426 White
<b>Product description</b>	 <p>VALUE folding shower seat and backrest create a simple solution for users who need to sit down while showering.</p> <p>Back rest for folding shower seat, back rest pale grey Incl. mounting kit for brick wall</p>
<b>Adjustability</b>	Fixed mounting
<b>Materials</b>	The wall profile is made of diecast aluminium, and the backrest is made of vacuum-moulded polyethylene (PE)
<b>Warranty</b>	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.



# PRESSALIT®

<b>Product</b>	VALUE backrest, fixed height
<b>Item number</b>	R1651000
<b>Cleaning</b>	<p>VALUE backrest.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use corrosive or abrasive cleaners as it may result in the surface being damaged.</p> <p>For surface disinfection, chlorine-based disinfectants may be used according to the product specification and safety precautions for these products.</p> <p>The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container.</p> <p><b>DISINFECTANTS:</b> Clinimax Difficil-S Erisan OXY+ Klorilli</p> <p><b>NON-DISINFECTANTS:</b> Suma Star free D1 Taski 100 Taski 200</p> <p>The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.</p>
<b>Declaration of conformity</b>	CE MDR This product is in compliance with all relevant requirements of: - Regulation (EU) 2017/745 on medical devices.

[pressalit.com](http://pressalit.com)

## PRESSALIT®