

## PRESSALIT.

| Product<br>Item number | CareSystem support arm, fixed height R3016000                                                                                                                      |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Colour                 | White                                                                                                                                                              |
| Product description    | Support arm without counter balance, 600 mm, fixed height, folds against the wall Incl. mounting kit for brick wall.  Discontinued as per 31 Dec 2015              |
|                        | Support arm without counter balance, 600 mm, fixed height, folds against the wall Incl. mounting kit for brick wall.  Discontinued as per 31 Dec 2015              |
| Adjustability          | Fixed mounting                                                                                                                                                     |
| Product capacity       | Max. vertical load: 135 kg.                                                                                                                                        |
| Materials              | Basic material: 20 micrometre anodized aluminium                                                                                                                   |
|                        | Surface material: Two-component plastic (PVC)                                                                                                                      |
|                        | End cap/cover cap: Arnite                                                                                                                                          |
| Ergonomy               | The oval shape gives an easy and comfortable surface to lean and push against.                                                                                     |
|                        | The end of the support arm is enlarged so that it is easier to pull forward from a sitting position and also provides a clear sense of where the support arm ends. |
|                        | Soft and slip-resistant material that has a warm feel.                                                                                                             |
| Features               | Folds against the wall.                                                                                                                                            |
|                        | The underside is shaped for the mounting of various accessories.                                                                                                   |



## PRESSALIT®

| Product<br>Item number | CareSystem support arm, fixed height R3016000                                                                                                                                                                                                                                                                            |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Cleaning               | Pressalit recommends that chalky deposits be removed from wall tracks and other aluminum parts using a solution of three parts of 30% acetic acid to seven parts tap water. When the aluminum parts have been washed down with this solution, the parts should be wiped over with a damp cloth wrung out in clear water. |

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com