




# PRESSALIT®

<b>Product</b> <b>Item number</b>	CareSystem shower seat 410, manually height adjustable R7143000
<b>Colour</b>          <b>Product description</b>	White             Shower chair 410 with aperture, "soft feel", height adjustable 250 mm, Excl. mounting kit. Discontinued as per 31/12-2015  Shower chair 410 with aperture, "soft feel", height adjustable 250 mm, Excl. mounting kit. Discontinued as per 31/12-2015
<b>Adjustability</b>  <b>Product capacity</b>	Height adjustable  Max. vertical load: 150 kg.
<b>Materials</b>          <b>Ergonomy</b>       <b>Features</b>	Basic material: 20 micrometre anodized aluminium  Seat: Polyurethane foam with marine ply core  Seat surface: Gum lac  Lever: Stainless steel w/nylon knob  Front plate/end covers, track: Polystyrene plastic (PS)  Soft and slip-resistant material with a warm feel.  Ergonomically designed seat with aperture that makes intimate hygiene easier.  The seat can fold up when not in use or to make more space.  Height adjustable 250 mm.



# PRESSALIT®

<b>Product</b>	CareSystem shower seat 410, manually height adjustable
<b>Item number</b>	R7143000
<b>Cleaning</b>	Pressalit recommends that chalky deposits be removed from wall tracks and other aluminum parts using a solution of three parts of 30% acetic acid to seven parts tap water. When the aluminum parts have been washed down with this solution, the parts should be wiped over with a damp cloth wrung out in clear water.

[pressalit.com](http://pressalit.com)

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

## PRESSALIT®