

PRESSALIT.

Product Item number	PLUS shower seat 450, manually height and manually sideways adjustabl R7454000
	VVS:777516200 NRF:6317129
	EAN: 5708590359885
	HMI: 124239
	FDA: INN/D135865
Colour	White
Product description	PLUS shower seats are offered in a number of colours. They are available as both fixed height, electrically adjustable and manually adjustable, with and without aperture, and with an option of back and armrests. They are suitable for wall-mounting or for fitting to Pressalit PLUS wall tracks.
	PLUS shower seat 450, with backrest. Manually height adjustable by 380 mm and sideways adjustable, without the use of tools. The seat is foldable. Materials: Comfortable polyurethane foam over a plywood core, fitted to an aluminic profile. For fitting to horizontal PLUS wall track (not included). <i>Replaces R7440</i>
Adjustability	Height and sideways adjustable
Adjustment height	380 mm
Product capacity	Maximum user mass: 200 kg. When the seat is loaded with 200 kg, each wall fastener is loaded with 23 kg.
	The seat has been tested statically with safety factor 1.5, which is a requirement according to the applied assistive products standards.
Materials	The shower seat is made of polyurethane (PU) foam with a plywood core and this is covered with an in-mould coated surface that is protected using a clear PU lacquer. The backrest is made of polyurethane (PU) foam protected using a clear PU lacque Internal components and fasteners are primarily made of stainless steel. The hinge barrel is capped with ABS.
Ergonomy	Ergonomically designed seat - contact surface are 450 x 437 mm. Soft and slip- resistant material that has a warm feel.



PRESSALIT.

Product	PLUS shower seat 450, manually height and manually sideways adjusta
Item number	R7454000
Features	The seat can fold up when not in use or to maximise floor space. The seat and backrest can be adjusted simultaneously thanks to a safe and easy click system.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.
Cleaning	
	PLUS shower seat.
	The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged.
	The following cleaning agents have been tested by Pressalit and may be used, a long as they are diluted according to the instructions on the container.
	DISINFECTANTS:
	Erisan OXY+
	Tvättsprit 70%
	CaviCide surface disinfectant decontaminant cleaner
	Rodalon Desinfektionsspray Terralin protect (SCHÜLKE+)
	Dismonzon Plus (BODE Chemie GmbH)
	Districtizen Plus (BODE Chemie Gribit)
	NON-DISINFECTANTS:
	Wetrok Caledor (Wetrok)
	Petrosol Citro (Büfa)
	Important: Polyurethane surfaces of PLUS shower seats and PLUS backrests contains an antibacterial surface, which will be compromised if cleaned with anyt other than soap and water.
	The product must be wiped down with a soft rag wrung out in clear water immedi after cleaning, as some materials may otherwise be discolored and stained, and surface coating may deteriorate.
	DECALCIFICATION OF ALUMINUM PARTS: Pressalit recommends that chalky deposits be removed from wall tracks and other
	aluminum parts using a solution of three parts 30% acetic acid to seven parts tap water. When the aluminum parts have been washed down with this solution, the must be wiped over with a damp cloth wrung out in clear water.
Declaration of conformity	CE MDR This product is in compliance with all relevant requirements of: - Regula (EU) 2017/745 on medical devices.



PRESSALIT.

Product Item number PLUS shower seat 450, manually height and manually sideways adjustable R7454000

pressalit.com

Pressalit A/S Pressalitvej 1 DK-8680 Ry Tel.: +45 8788 8788 Fax: +45 8788 8789 E-mail: pressalit@pressalit.com

