

PRESSALIT.

Product	VALUE II fold down grab bar, 33.5" R1181000
item number	
	VVS:615892100
	NRF:6191059
	EAN:5708590373010
	HMI: 138203
Color	White
Product description	VALUE and VALUE II fold down grab bars with the ALBA wing mounting profile are specifically designed for use with SELECT TL3 and TL4 toilet lifters. The height adjustment of the toilet lifter is done only with the control panel on the fold down grab bar.
	VALUE II fold down grab bar, 33.5", foldable, height adjustable 5.9" without the use of tools (right hand operated). Provides firm support when sitting, standing, or transferring to and from a toilet. Materials: Rilsan-coated steel and diecast aluminum. For fitting to TL3/TL4 electrical/manual toilet lifter, ALBA wing mounting profile included. For mounting on user's left side.
Adjustability	Height adjustable
Adjustment height	5,91 inch
Product capacity	Maximum load: 265 lbs
Materials	The fold down grab bar is manufactured from steel with a Rilsan coating, which is comfortable to the touch, easy to clean, and protects the metal against corrosion. The ALBA wing mounting profile is made from diecast aluminium.
Ergonomy	VALUE II support arm with ALBA wing, 850 mm, white. Height adjustable, right hand operated, for fitting to TL3/TL4 (left)
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.



PRESSALIT.

Product	VALUE II fold down grab bar, 33.5"
Item number	R1181000
Cleaning	
J	This product should be cleaned with hot water or a mild cleaning solution. Do not use corrosive or abrasive cleaners as it may result in the surface being damaged.
	For surface disinfection, chlorine-based disinfectants may be used according to the product specification and safety precautions for these products.
	The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container.
	DISINFECTANTS:
	Clinimax Difficil-S
	Erisan OXY+ Sprit (husholdningssprit)
	Suma Bac D10
	Virkon
	Rodalon Desinfektionsspray
	Lysoformin Spezial CaviCide
	Sufanios Premium
	Tvättsprit 70%
	WipeClean UniversalCleaning
	Chrolorine Disinfection Wipes
	Klorilli (Kiilto Clean)
	NON-DISINFECTANTS:
	Suma Star free D1
	Taski 100
	Taski 200
	Milizid Citro
	The product must be wiped down with a soft rag wrung out in clear water immediat after cleaning, as some materials may otherwise be discolored and stained, and th surface coating may deteriorate.
Declaration of conformity	This product is in compliance with all relevant requirements of regulation (EU) 2017/745 on medical devices.



PRESSALIT®

Product VALUE II fold down grab bar, 33.5"

Item number R1181000

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com