

PRESSALIT.

Product	Mattress for MCT 2 and MCT 3
Item number	R842717378
	VVS:777542718
	EAN:5708590373881
	HMI: 128440
Colour	Navy blue
Product description	Mattress for MCT 2 and MCT 3 change trolleys, 695 x 1680 mm. Can also be used with MSCT 1
	Mattress for MCT 2 and MCT 3 change trolleys, 695 x 1680 mm.
Materials	The mattress cover is manufactured from 100% polyester with polyurethane, approved in accordance with Oeko-Tex 100, Class I. It is fire retardant and impermeable to dust mites. The foam filling is flexible polyester-based polyurethan foam, which is also fire retardant and approved in accordance with Oeko-tex 100, Class I.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of three years from the invoicing date.



PRESSALIT®

Product	Mattress for MCT 2 and MCT 3
ltem number	R842717378
Cleaning	Mattress for CT and MCT. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged. The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container. DISINFECTANTS: Sprit (husholdningssprit) Tvättsprit 70% Safe Zone Plus Selgiene Extreme Accel PREVention Wipes Lingettes The product must be wiped down with a soft rag wrung out in clear water immedia after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com