

| Product | SELECT TL3 toilet lifter, manual, for floor outlet | |
|---------------------|---|--|
| Item number | R9013000 | |
| | VVS:617989200 | |
| | NRF:6191177 | |
| | EAN:5708590372754 | |
| | HMI: 138231 | |
| Color | White | |
| Product description | SELECT TL3 and TL4 toilet lifters are developed with a clear focus on hygiene. Please contact Pressalit for more information about compatible toilet bowls for Pressalit SELECT toilet lifters | |
| | SELECT TL3 toilet lifter, for floor outlet. Incl. cUPC-approved cistern with dual flushing system preset at 0.8/1.6 gallons (can be changed to 0.5/1.1 gallons). Manually height adjustable 8" with crank handle. Toilet pan, toilet seat and fold down grab bars to be ordered separately. Materials: Stainless steel and PMMA/ABS. Excl. wall mounting kit. | |
| Adjustability | Height adjustable | |
| Adjustment height | 7,88 inch | |
| Product capacity | Maximum user weight: 440 lbs. Maximum load: 518 lbs. (user weight + accessories). When the toilet lifter is loaded with 518 lbs, each wall fastener is loaded with 143 lbs. The toilet lifter has been tested statically with 1320 lbs, which is a requirement according to the applied assistive products standards. Please note that most toilets compatible with this lifter have been tested up to 880 lbs. | |
| Materials | The toilet lifter is constructed of stainless steel, the cover fairing is made of PMMA/ABS plastics. | |
| Warranty | Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date. | |



| Product Item number | SELECT TL3 toilet lifter, manual, for floor outlet R9013000 |
|------------------------|--|
| Cleaning | This product should be cleaned with hot water or a mild cleaning solution. Do not use corrosive or abrasive cleaners as it may result in the surface being damaged. For surface disinfection, chlorine-based disinfectants may be used according to the product specification and safety precautions for these products. The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container. DISINFECTANTS: CaviCide The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate. |



| Product | | | | |
|---------|--------|--|--|--|
| Item | number | | | |

SELECT TL3 toilet lifter, manual, for floor outlet

R9013000

Mounting in different wall types

| Concrete >B25 V8605 - Mounting kit, 4 lag bolts Ø10 x 70, plugs SX12, with grey covers. To be ordered separately. | |
|---|--|
| Solid Brick >= Mz12 V8605 - Mounting kit, 4 lag bolts Ø10 x 70, plugs SX12, with grey covers. To be ordered separately. | |
| Solid Sand Lime Brick >= HIz12 V8605 - Mounting kit, 4 lag bolts Ø10 x 70, plugs SX12, with grey covers. To be ordered separately. | |
| Aerated Concrete P4 V8655 - Mounting kit, 4 lag bolts Ø10 x 70, plugs GB14, with grey covers. To be ordered separately. | |
| Drywall (2 x 12,5 mm) water-resistant, re-enforced with wood V8605 - Mounting kit, 4 lag bolts Ø10 x 70, plugs SX12, with white covers. NB! Plugs are not needed. To be ordered separately. | |



| Product Item number | SELECT TL3 toilet lifter, manual, for floor outlet R9013000 | |
|------------------------|---|---|
| | Drywall (2 x 12,5 mm) water-resistant, re-enforced with 22 mm plywood Wall fasteners (4 pcs) must be bought separately, for drywall it is possible to use e.g. SnapTight or WingIt anchors or similar. | |
| | Drywall (2 x 12,5 mm) water-resistant, re-enforced with steel plates (.08 inch) Wall fasteners (4 pcs) must be bought separately, for drywall it is possible to use e.g. SnapTight or WingIt anchors or similar. | |
| eransinal (| Drywall (2 x 12,5 mm) water-resistant | _ |
| | Mounting in this wall type is not recommended - risk of fatal injury to persons and walls. | |

pressalit.com