




# PRESSALIT®

<b>Product</b> <b>Item number</b>	Wall-hung toilet 700 mm R2090000
<b>Colour</b>	VVS:613231000 NRF:6191648 EAN:5708590207971 HMI: 35272 NBS:N13:10,11,300 White
<b>Product description</b>	 <p><b>Discontinued per 31 Dec 2018 and replaced by R2099</b></p> <p>Toilet, wallhung, distance from wall to the front: 700 mm Excl. mounting kit</p> <p><b>Discontinued per 31 Dec 2018 and replaced by R2099</b></p> <p>Toilet, wallhung, distance from wall to the front: 700 mm. Excl. mounting kit</p>
<b>Product capacity</b>	Max load: 400 kg.
<b>Materials</b>	Ceramic
<b>Ergonomy</b>	700 mm is recommended for wheelchair users as it gives more space for a carer. Wallhung toilet makes mounting possible at individual height. Leaves space for foot support, the legs of a mobile hoist, and the carer's feet.
<b>Warranty</b>	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.



# PRESSALIT®

<b>Product</b>	Wall-hung toilet 700 mm
<b>Item number</b>	R2090000

---

[pressalit.com](http://pressalit.com)

---

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

# PRESSALIT®