




# PRESSALIT®

<b>Product</b>	Wall-hung toilet 540 mm
<b>Item number</b>	R2091000
<b>Colour</b>	VVS:613230000 NRF:6190522 EAN:5708590213897 HMI: 35273 NBS:N13:10,11,300 White
<b>Product description</b>	 <p>Toilet, wall-hung, distance from wall to the front: 540 mm Excl. mounting kit</p> <p>Toilet, wall-hung, distance from wall to the front: 540 mm. Excl. mounting kit</p>
<b>Product capacity</b>	Max load: 400 kg.
<b>Materials</b>	Ceramic
<b>Ergonomy</b>	Wall-hung toilet makes mounting possible at individual height. Leaves space for foot support, the legs of a mobile hoist, and the carer's feet.
<b>Warranty</b>	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.

[pressalit.com](http://pressalit.com)

# PRESSALIT®