




PRESSALIT®

Product Item number	Freestanding column 945 mm, for height adjustable PLUS support arm R9431
Product description	VVS:615856990 NRF:6190668 EAN:5708590320434 NBS:N13:436  Free-standing counn, grey, for height adjustable PLUS support arm Free-standing counn, grey, for height adjustable PLUS support arm
Adjustment height	240 mm
Materials	Basic material: Steel Surface material: Rilsan
Features	Used with freestanding toilets and where it is not possible to wallmount a support arm. For adjustable height support arm.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.



PRESSALIT®

Product	Freestanding column 945 mm, for height adjustable PLUS support arm
Item number	R9431

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®

8. juli 2026