




Product Item number	Pressalit Public Universal Disinfection dispenser with sensor, 950 ml, brushed steel PZ30732
Product description	<p>VVS: 774566615 NRF: 1438660 RSK: 8920592 EAN: 5708590396361</p>  <p>Disinfection dispenser with stainless steel sensor for wall mounting. One-handed disinfection pump with 950 ml tank. Designed for standard hand sanitizers. Lockable solution. Including mounting kit.</p> <p>Inspection opening indicates filling level.</p>
Materials	European stainless steel (AISI 316L)
Warranty Cleaning WEEE	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Please enquire of your local Pressalit subsidiary or distributor (in countries where there is no subsidiary) that will arrange for recycling of electronic waste in accordance with the WEEE directive.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>



PRESSALIT®

Product	Pressalit Public Universal Disinfection dispenser with sensor, 950 ml, brushed steel
Item number	PZ30732

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®