




PRESSALIT®

Institutional hinge, universal

Item number

BV5999
EAN:5708590267265



*55 mm

Product capacity

Weight

0,257 kg

Gross weight

0,257 kg

Weight per crate

0 kg

Qty per crate

0

Qty per pallet

0

The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.

Subject to changes in the product range, technical changes and printing errors.

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®