




# PRESSALIT®

|                            |                                                                                                                                                                                                 |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product</b>             | Mixer tap                                                                                                                                                                                       |
| <b>Item number</b>         | RT633                                                                                                                                                                                           |
| <b>Product description</b> | <p>RSK:8064034<br/>EAN:5708590396309</p>  <p>Mixer tap for daycare changing tables and basin vanity unit</p> |
| <b>Materials</b>           | Brass with chrome plating                                                                                                                                                                       |
| <b>Warranty</b>            | Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.                                                                   |
| <b>Cleaning</b>            | The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.                      |

[pressalit.com](http://pressalit.com)