



PRESSALIT®

Product	Mounting kit V8657 (9 pcs), for aerated concrete
Item number	V8657
Product description	<p>EAN:5708590351377</p> <p>Mounting kit, 9 lag bolts Ø10 x 70, plugs GB14, for aerated concrete.</p> <p>Used for: INDIVO motorized lifts for kitchen counter-top, wall-mounted RK1810 RK1811 RK1812 RK1813 RK1814 RK1850 RK1851 RK1852 RK1853 RK1854</p> <p>Mounting kit, 9 lag bolts Ø10 x 70, plugs GB14, for aerated concrete.</p>
Materials	Stainless steel and plastics



PRESSALIT®

Product	Mounting kit V8657 (9 pcs), for aerated concrete
Item number	V8657

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®