



PRESSALIT®

Product	Mounting kit V8654 (4 pcs), for aerated concrete
Item number	V8654
	EAN:5708590351407
Product description	Mounting kit, 4 countersunk bolts M8 x 35, plugs FTP M 8, for aerated concrete. Used for: PLUS fold down grab bars: R360070, R360085 PLUS folding seat: R5510, R5504 PLUS shower seats: R7304, R7404 Mounting kit, 4 countersunk bolts M8 x 35, plugs FTP M 8, for aerated concrete
Materials	Stainless steel

pressalit.com