



# PRESSALIT®

## Pressalit Trend incl. Pressalit Sway D2 994

Wall hung toilet with circular flush, incl. toilet seat  
Toilet seat with soft close, lift-off and antibacterial effect, incl. hinge in stainless steel

### Item number



CP994011CF  
VVS: 613326500  
NRF: 1462969  
RSK: 7822116  
EAN: 5708590396903

The Pressalit Trend Combi Pack with Sway D2 seat offers a perfect union of minimalist elegance and advanced functionality. The rimless, wall-hung toilet features the hygienic and efficient Circular Flush system for a quiet, splash-free flush. The award-winning Sway D2 seat is slim and stylish, complementing the bowl's modern design while providing excellent comfort. Soft Close ensures a gentle, noiseless close, and Lift-Off makes cleaning simple. An antibacterial additive helps reduce bacterial growth for improved hygiene. This combination is ideal for contemporary bathrooms where design and practicality are equally important, creating a refined, easy-to-maintain solution for everyday use.




### Colour

White



\*20 mm / \*\* 40 mm





---

**Pressalit Trend incl. Pressalit Sway D2 994**

**Wall hung toilet with circular flush, incl. toilet seat**  
**seal with soft close, lift-off and antibacterial effect, incl. hinge in stainless steel**

---

**Materials****TOILET SEAT:**

The material is colour ingrained duroplast (UF A 10 = urea formaldehyde) that contains no environmentally hazardous substances. UF plastic comprises 67% urea formaldehyde resin, 28% cellulose, as well as 5% minerals, pigments, lubricants and moisture content. The raw material from the supplier contains up to 0.03% free formaldehyde and even less following curing.

**ANTIBACTERIAL:**

This toilet seat features an antibacterial additive with biocidal silver chloride, which helps prevent bacterial growth.

**BUFFERS:**

EVA (copolymer made of ethylene and vinyl acetate).

**DAMPER KIT (only relevant for toilet seats with soft close):**

Hydraulic damper with stainless steel damper casing and thermoplastic parts.

**MOUNT/HINGES:**

Stainless steel model W. no. 1.4301/M and thermoplastic.

**Product capacity**

Load - ring seat 240 kg

**Weight**

32,035 kg

**Carton size**

540x410x510 mm

**Gross weight**

34,704 kg

**Weight per crate**

34,704 kg

**Qty per crate**

1

**Crates per pallet**

5

**Qty per pallet**

5

**Warranty**

Pressalit provides a 10-year warranty on toilet seats and a 5-year warranty on toilet bowls against presentation of a valid receipt. In the event of a complaint, please contact the retailer or store where the product was purchased. The warranty covers defective materials and/or manufacturing defects and/or functional limitations that may occur following proper use and cleaning of the product. Pressalit reserves the right to examine the product at our factory. In the case of a valid complaint, Pressalit will either:  
- supply spare parts, or  
- replace the product with an identical or a similar model.  
If the color of the toilet seat is no longer available, Pressalit reserves the right to supply a similar model in white. Upon replacement or delivery of spare parts, the warranty period will not be extended nor renewed. If replacement is not possible, Pressalit is entitled instead to refund the purchase price.



# PRESSALIT®

---

**Pressalit Trend incl. Pressalit Sway D2 994**

Wall hung toilet with circular flush, incl. toilet seat  
Toilet seat with soft close, lift-off and antibacterial effect, incl. hinge in stainless steel

---

Subject to changes in the product range, technical changes and printing errors.

---

[pressalit.com](http://pressalit.com)

---

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

# PRESSALIT®