




PRESSALIT®

Product Item number	Pressalit Public Classic Disinfection dispenser, 950 ml, polished steel PY10730
Product description	VVS: 774565540 NRF: 1438620 RSK: 8920384 EAN: 5708590392288  Stainless steel disinfection dispenser for wall mounting. One-hand operated disinfection pump with 950 ml tank. Designed for standard hand sanitizers. Lockable solution. Including mounting kit. Inspection opening indicates filling level.
Materials	European stainless steel (AISI 316L)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com