




PRESSALIT®

Product Item number	Pressalit Public Universal Disinfection dispenser, 950 ml, brushed steel PZ30730
Product description	VVS: 774566515 NRF: 1438659 RSK: 8920591 EAN: 5708590396354  Stainless steel disinfection dispenser for wall mounting. One-hand operated disinfection pump with 950 ml tank. Designed for standard hand sanitizers. Lockable solution. Including mounting kit. Inspection opening indicates filling level.
Materials	European stainless steel (AISI 316L)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com