



PRESSALIT®

Product Item number	Pressalit Choice Towel holder, single, 60 x 8,5 cm, brushed brass Q70500
Product description	VVS:775711105 NRF:7001082 RSK:8755026 EAN:5708590372235 DB: 2173272 Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.
Materials	Stainless steel AISI 304, mounting knot in brass
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com