




PRESSALIT®

Product Item number	Pressalit Public Premium Soap dispenser, 200 ml, brushed steel PX30710
Product description	VVS:774684215 NRF:1438645 RSK:8920370 EAN:5708590392172  Stainless steel soap dispenser for wall mounting. One-hand operated and with drip-free soap pump with a 200 ml container. Designed for regular liquid soap or hand lotion. Incl. mounting kit.
Materials	European stainless steel (AISI 316L)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com