




Product Item number	TMA 1 arm support for mounting on toilet R1171000
Colour Product description	<p>VVS:615723100 NRF:6191651 EAN:5708590211718 HMI: 42230 NBS:N13:436 DB: 2076624</p> <p>White</p>  <p>TMA arm supports comes in various widths, and models with or without toilet seat.</p> <p>TMA 1 arm supports for mounting on toilet, incl. BL1999 hinge. Dark grey handles. BL1 hinge, is only recommended with toilet seats R26, R27, R28, R29, R37 or R38.</p>
Adjustability	Fixed mounting
Materials Ergonomy Features	<p>Powder-coated stainless steel. Soft polyurethane. <i>Flammability-tested (ISO 8191-1). Test report available.</i></p> <p>The folding arms give a good support and the moulded handles are comfortable for the hands and arms.</p> <p>The arm supports can be folded individually.</p>
Warranty Cleaning Declaration of conformity	<p>Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.</p> <p>GENERAL The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged.</p> <p>The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.</p> <p>CE UK MDR This product is in compliance with all relevant requirements of: - Regulation (EU) 2017/745 on medical devices. For the UK: - Medical Device Regulation (MDR) 2002 on medical devices.</p>



PRESSALIT®

Product	TMA 1 arm support for mounting on toilet
Item number	R1171000

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®