




PRESSALIT®

Product Item number	Pressalit Public Classic Soap dispenser, 1100 ml, polished steel PY10710
Product description	VVS: 774565140 NRF: 1438644 RSK: 8920381 EAN: 5708590392257  Stainless steel soap dispenser for wall mounting. One-handed, drip-free soap pump with 1100 ml container. Designed for regular liquid soap or hand cream. Lockable solution. Incl. mounting kit.
Materials	European stainless steel (AISI 316L)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com