



PRESSALIT®

Product Item number	Changing table mattress 950 mm R8703378
Colour	VVS:777630808 RSK:8565132 EAN:5708590366081 HMI: 125808 Navy blue 
Product description	Mattresses for daycare changing tables are available in several colors and in lengths that are adapted to the different changing table lengths. Accessories for changing tables. Mattress 760 mm x 950 mm x 30 mm, for 1400 mm changing table with sanitary appliances. Phthalate-free
Materials	Cover: Techmadin, phthalate-free Padding: Foam
Cleaning	Upholstery (Cover): Techmadin, phthalate-free Daily cleaning with soapy water and subsequent wiping over with a damp cloth, and wash regularly at 60-90°, without fabric softener. Shrinkage after washing may occur. Do not turn inside out before washing. Can be tumble-dried at max. 60°. Can be disinfected with antiseptic wipes. Alcohol % < 70. Do not use disinfectants intended for hard surfaces. Padding: (foam) The foam contains only DINP softeners which are approved by the EU for the use in articles for children. Do not wash the foam. Do not use the mattress without the cover for reasons of hygiene.



PRESSALIT®

Product	Changing table mattress 950 mm
Item number	R8703378

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®

24. september 2026