




PRESSALIT®

Product Item number	Mounting kit V8658 (6 pcs),for aerated concrete V8658
Product description	EAN:5708590351360  Mounting kit, 6 coach screws Ø10 x 70, plugs GB14, for aerated concrete. Used for: SCT 1000 shower change tables: R8402, R8403 SCT 2000 shower change tables: R8478, R8582418 R8572418(AU) SCT 3000 shower change tables: R8582313, R8582318 R8572313(AU), R8572318(AU) Changing tables, low start, wallhung: R8661, R8662, R8663, R8671, R8672, R8673 Changing tables, wallhung: R8732, R8734, R8733, R8735, R8711, R8712, R8713, R8721, R8722, R8723 INDIVO Lift for wall cupboard: RK1220, RK1221 (AU), RK1420, RK1421 MATRIX basin unit, height and sideways adjustable: R482519(UK), R482719(UK), R484519(UK), R484719(UK), R485519(UK), R485719(UK), R487519(UK) R482019, R482219, R484019, R484219, R485019, R485219, R487019 Mounting kit, 6 coach screws Ø10 x 70, plugs GB14, for aerated concrete
Materials	Stainless steel and plastic



PRESSALIT®

Product	Mounting kit V8658 (6 pcs),for aerated concrete
Item number	V8658

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®

8. juli 2026