



Product	Mounting kit V8658 (6 pcs), for aerated concrete
Item number	V8658
Product description	<p>EAN:5708590351360</p> <p>Mounting kit, 6 lag bolts Ø10 x 70, plugs GB14, for aerated concrete.</p> <p>Used for: SCT 1000 shower change tables: R8402, R8403 SCT 2000 shower change tables: R8478, R8592418 SCT 3000 shower change tables: R8592313, R8592318 Changing tables, wall-mounted, low start: R8661, R8662, R8663, R8671, R8672, R8673 Changing tables, wall-mounted: R8732, R8734, R8733, R8735, R8711, R8712, R8713, R8721, R8722, R8723 INDIVO Lift for wall cabinet: RK1820, RK1821 MATRIX sink brackets, vertical and horizontal adjustable: R482619, 482819</p> <p>Mounting kit, 6 lag bolts Ø10 x 70, plugs GB14, for aerated concrete</p>
Materials	Stainless steel and plastics



PRESSALIT®

Product	Mounting kit V8658 (6 pcs), for aerated concrete
Item number	V8658

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®