




PRESSALIT®

Product Item number	Pressalit Public Premium Toilet paper holder for 1 roll of toilet paper + 1 spare roll, polished steel PX10210
Product description	<p>VVS: 776130240 NRF: 1438648 RSK: 8920358 EAN: 5708590392073</p>  <p>Toilet paper holder in steel housing with space for spare roll. For wall mounting. The holder is lockable and includes mounting kit.</p>
Materials	European stainless steel (AISI 304)
Warranty Cleaning	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>

pressalit.com

PRESSALIT®