




PRESSALIT®

Product Item number	TMA 1 arm support for mounting on toilet R1170000
Colour	VVS:615723200 NRF:6191649 RSK:7828180 EAN:5708590211701 HMI: 35270 NBS:N13:436 FDA:ILS/D135866 White
Product description	 TMA arm supports comes in various widths, and models with or without toilet seat. TMA 1 arm supports for mounting on toilet, incl. R37 Dania toilet seat. With BL1999 hinge. Dark grey handles.
Adjustability	Fixed mounting
Materials	Powder-coated stainless steel. Soft polyurethane. <i>Flammability-tested (ISO 8191-1). Test report available.</i> Toilet seat and cover in thermoset plastic (urea).
Ergonomy	The folding arms give a good support and the moulded handles are comfortable for the hands and arms.
Features	The arm supports can be folded individually.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.
Cleaning	The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged. The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.



PRESSALIT®

Product	TMA 1 arm support for mounting on toilet
Item number	R1170000

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®

24. september 2026