




PRESSALIT®

Product Item number	Mounting kit V8652 (2 pcs), for aerated concrete V8652
Product description	<p>EAN: 5708590366005</p>  <p>Mounting kit, 2 coach screws Ø10 x 100, plugs GB14, for aerated concrete.</p> <p>Used for: MATRIX SMALL wash basins: R2020/R2022, R2021/R2023 MATRIX CURVE wash basins: R2050, R2051, R2052, R2053</p> <p>Mounting kit, 2 coach screws Ø10 x 100, plugs GB14, for aerated concrete</p>
Materials	Stainless steel and plastic

pressalit.com