



PRESSALIT®

Product	Pressalit Public Classic Disinfection dispenser with sensor, 700 ml, brushed steel
Item number	PY30725
Product description	<p>VVS: 774565615 NRF: 1438614 RSK: 8920390 EAN: 5708590392745</p> <p>Disinfection dispenser with stainless steel sensor for wall mounting. One-handed disinfection pump with 700 ml tank. Designed for standard hand sanitizers. Lockable solution. Including mounting kit.</p> <p>Inspection opening indicates filling level.</p>
Materials	European stainless steel (AISI 316L)
Warranty	A 2-year warranty is provided for material and manufacturing defects.
Cleaning	This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.
WEEE	<p>Please enquire of your local Pressalit subsidiary or distributor (in countries where there is no subsidiary) that will arrange for recycling of electronic waste in accordance with the WEEE directive.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>

pressalit.com