



PRESSALIT®

<b>Product</b>	Wall-mounted toilet 27 4/8"
<b>Item number</b>	R2099000
<b>Color</b>	VVS:612928000 NRF:6191026 RSK:7822154 EAN:5708590345284 HMI: 105700 White
<b>Product description</b>	Toilet, wall-mounted, distance from wall to the front: 27 4/8" Excl. mounting kit. For mounting directly on wall, a mounting kit from V&B must be ordered from Pressalit at an additional cost. For mounting on built-in cistern or Pressalit toilet lifter, mounting kit is included in the delivery.  Toilet, wall-mounted, distance from wall to the front: 27 4/8". Excl. mounting kit.
<b>Product capacity</b>	Max load: 880 lbs.
<b>Materials</b>	Sanitary porcelain
<b>Ergonomy</b>	27 4/8" is recommended for wheelchair users as it provides sufficient space for a caregiver. Wall mounted toilet makes mounting possible at individual height. Leaves space for foot support, the legs of a floor lift, and the caregiver's feet.
<b>Warranty</b>	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.



# PRESSALIT®

<b>Product</b>	Wall-mounted toilet 27 4/8"
<b>Item number</b>	R2099000

---

[pressalit.com](http://pressalit.com)

---

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

PRESSALIT®