




# PRESSALIT®

|                                      |                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|--------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product</b><br><b>Item number</b> | Pressalit Public Classic Disinfection dispenser with sensor, 700 ml, matt black<br><br>PY40725                                                                                                                                                                                                                                                                                                                                      |
| <b>Product description</b>           | <p>VVS: 774565631<br/>NRF: 1438615<br/>RSK: 8920399<br/>EAN: 5708590392769</p>  <p>Disinfection dispenser with stainless steel sensor for wall mounting. One-handed disinfection pump with 700 ml tank. Designed for standard hand sanitizers. Lockable solution. Including mounting kit.</p> <p>Inspection opening indicates filling level.</p> |
| <b>Materials</b>                     | European stainless steel (AISI 316L)                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Warranty</b><br><b>Cleaning</b>   | <p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>                                                                               |

[pressalit.com](http://pressalit.com)