



PRESSALIT®

Product Item number	Pressalit Public Classic Waste bin 60 l, polished steel PY10920
Product description	VVS:777744640 NRF:1438613 RSK:8920386 EAN:5708590392301 Stainless steel waste bin for floor mounting or concealed wall mounting. Front panel with a stainless steel plate with rounded corners. Equipped with removable, washable polypropylene fabric bag. Capacity approx. 60 liters. Comes with mounting kit. 2.5 mm and 4 mm thick grating and mounting holes placed in the second row for concealed wall mounting.
Materials	European stainless steel (AISI 304)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com