



PRESSALIT®

Product Item number	Pressalit Choice Shelf, brushed brass Q70410
Product description	VVS:772471305 NRF:7001576 RSK:8890002 EAN:5708590381879 DB: 2425189  Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.
Materials	Stainless steel AISI 304
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com

PRESSALIT®