



PRESSALIT®

Splashguard, transparent mat, for 50 mm raised toilet seats

Item number	T9516 VVS: 616810516 NRF: 6190363 EAN: 5708590183909	
Description	Splashguard, transparent mat, for 50 mm raised toilet seats. Splashguard, transparent mat, for 50 mm raised toilet seats.	
Materials	Polypropylene (PP)	
Product capacity		
Weight	0,5 kg	
Gross weight	0,5 kg	
Weight per crate	0 kg	
Qty per crate	0	
Qty per pallet	0	

pressalit.com
