



United
Kingdom
Conformity
Assessment

United Kingdom RP:
Pressalit Ltd.
100, Longwater Avenue
Green Park
Reading RG2 6GP
United Kingdom
T: 01180923 2696
uk@pressalit.com

The undersigned hereby declares that the below-mentioned products comply with the relevant requirements of the following regulation/directives:
- Medical Device Regulation 2002
- The Supply of Machinery (Safety) Regulations 2008.

The products are manufactured in accordance with the following standards:	DS/ISO 17966:2016
---	-------------------

Model	R8012000	GMDN	63053	EAN	5708590345291
Description					
SELECT TL1 toilet bracket, for floor outlet, white cover plate. Manually height adjustable with crank handle					

Model	R8012112	GMDN	63053	EAN	5708590345307
Description					
SELECT TL1 toilet bracket, for floor outlet, anthracite cover plate. Manually height adjustable with crank handle					

Model	R8022000	GMDN	63053	EAN	5708590345321
Description					
SELECT TL1 toilet bracket, for wall outlet, white cover plate. Manually height adjustable with crank handle					

Model	R8022112	GMDN	63053	EAN	5708590345338
Description					
SELECT TL1 toilet bracket, for wall outlet, anthracite cover plate. Manually height adjustable with crank handle					

Ry, 5/17/2022  Trine Danielsen, Regulatory Consultant, Pressalit A/S

Signed by the person, who is authorized to compile the technical file, and empowered by the manufacturer to draw up this declaration.



Manufacturer:
Pressalit A/S - Pressalitvej 1 – DK-8680 Ry - T: +45 8788 8788 - pressalit@pressalit.com