



# PRESSALIT®

## Pressalit Trend incl. Pressalit Projecta D Solid Pro 1008

Wall hung toilet with circular flush, incl. toilet seat  
Toilet seat with soft close and antibacterial effect, incl. Solid Pro hinge

### Item number




CP1008011CF  
VVS: 613326300  
NRF: 1462950  
RSK: 7822114  
EAN: 5708590396880

The Pressalit Trend Combi Pack with Projecta D Solid Pro seat is designed for maximum durability and stability. Featuring the innovative rimless Trend bowl with Circular Flush for efficient, splash-free performance, it's the ideal solution for high-demand environments. The Solid Pro seat's robust D-shaped design and reinforced hinges withstand heavy, frequent use and tested with up to 60,000 sideways load cycles. An antibacterial additive helps maintain a hygienic surface and offers uncompromising quality. This combination is perfect for public buildings, institutions, or other settings where strength, stability, and hygiene are paramount, offering a reliable and long-lasting toilet solution.



### Colour

White





---

**Pressalit Trend incl. Pressalit Projecta D Solid Pro 1008**

**Wall hung toilet with circular flush, incl. toilet seat**  
**seal with soft close and antibacterial effect, incl. Solid Pro hinge**

---

**Materials**

**BOWL:**  
Porcelain

**TOILET SEAT:**

The material is colour ingrained duroplast (UF A 10 = urea formaldehyde) that contains no environmentally hazardous substances. UF plastic comprises approx. 70% urea formaldehyde resin, 30% cellulose, as well as 5% minerals, pigments, lubricants and moisture content. Free formaldehyde is released during the production process. Pressalit toilet seats have been measured after curing (< 0.01 mg/m<sup>3</sup>) to emit significantly below the permitted EU limit of 0.08 mg/m<sup>3</sup>.

**ANTIBACTERIAL:**

This toilet seat features an antibacterial additive with biocidal silver phosphate glass, which helps prevent bacterial growth.

**BUFFERS:**

EVA (copolymer made of ethylene and vinyl acetate).

**DAMPER KIT (only relevant for toilet seats with soft close):**

Hydraulic damper with stainless steel damper casing and thermoplastic parts.

**MOUNT/HINGES:**

Stainless steel. Thermoplastic semi-crystalline polyamides with partially aromatic content. Aluminium profile EN AW-6063-TS.

**Product capacity**

Load - ring seat 240 kg

**Weight**

31,508 kg

**Carton size**

540x410x510 mm

**Gross weight**

34,223 kg

**Weight per crate**

34,223 kg

**Qty per crate**

1

**Crates per pallet**

5

**Qty per pallet**

5



# PRESSALIT®

---

**Pressalit Trend incl. Pressalit Projecta D Solid Pro 1008**

**Wall hung toilet with circular flush, incl. toilet seat  
Toilet seat with soft close and antibacterial effect, incl. Solid Pro hinge**

---

## Warranty

Pressalit provides a 10-year warranty on toilet seats and a 5-year warranty on toilet bowls against presentation of a valid receipt. In the event of a complaint, please contact the retailer or store where the product was purchased. The warranty covers defective materials and/or manufacturing defects and/or functional limitations that may occur following proper use and cleaning of the product. Pressalit reserves the right to examine the product at our factory. In the case of a valid complaint, Pressalit will either:  
- supply spare parts, or  
- replace the product with an identical or a similar model.  
If the color of the toilet seat is no longer available, Pressalit reserves the right to supply a similar model in white. Upon replacement or delivery of spare parts, the warranty period will not be extended nor renewed. If replacement is not possible, Pressalit is entitled instead to refund the purchase price.

Subject to changes in the product range, technical changes and printing errors.

---

[pressalit.com](http://pressalit.com)

---

PRESSALIT®