



# PRESSALIT®

## Raised hinge, 100 mm


**Item number**

T9801

**Description**

EAN:5708590183916

Raised hinge, 100 mm

**Product capacity****Weight**

0,46 kg

**Gross weight**

0,46 kg

**Weight per crate**

0 kg

**Qty per crate**

0

**Qty per pallet**

0

This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.

[pressalit.com](http://pressalit.com)

# PRESSALIT®