



United Kingdom RP: Pressalit Ltd. 100, Longwater Avenue Green Park Reading RG2 6GP United Kingdom T: 01180923 2696 uk@pressalit.com

The undersigned hereby declares that the below-mentioned products comply with the relevant requirements of the following regulation/directives:

- Medical Device Regulation 2002

- The Supply of Machinery (Safety) Regulations 2008.

	DS/ISO 17966:2016
The products are manufactured in accordance with the following standards:	

Model	R9014000	GMDN	63053	EAN	5708590372914	
Description						
SELECT TL4 toilet bracket, for floor outlet, prepared for shower toilet seat, white single molded cover. Manually height						
adjustable with crank handle						

Model	R9014112	GMDN	63053	EAN	5708590372907	
Description						
SELECT TL4 toilet bracket, for floor outlet, prepared for shower toilet seat, anthracite single molded cover. Manually						
height adjustable with crank handle						

Model	R9024000	GMDN	63053	EAN	5708590372990
Description					
SELECT TL4 toilet bracket, for wall outlet, prepared for shower toilet seat, white single molded cover. Manually height					
adjustable with crank handle					

Model	R9024112	GMDN	63053	EAN	5708590372983
Description					
SELECT TL4 toilet bracket, for wall outlet, prepared for shower toilet seat, anthracite single molded cover. Manually					
height adjustable with crank handle					

Ry, 5/17/2022

Trine Danielsen, Regulatory Consultant, Pressalit A/S

Signed by the person, who is authorized to compile the technical file, and empowered by the manufacturer to draw up this declaration.



## Manufacturer:

Pressalit A/S - Pressalitvej 1 – DK-8680 Ry - T: +45 8788 8788 - pressalit@pressalit.com