



PRESSALIT®

Product Item number	Pressalit Public Premium Hygienic waste bin 4 l, brushed steel PX30925
Product description	VVS:777743415 NRF:1438634 RSK:8920372 EAN:5708590392196 Wall-mounted sanitary bin with bag dispenser integrated on the side, removable 4 liter sanitary bin and hinged lid. Includes mounting kit and a box with plastic sanitary bags.
Materials	European stainless steel (AISI 304)
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com